

Product datasheet for **MR204834L3V**

Tal1 (NM_011527) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tal1 (NM_011527) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tal1
Synonyms:	bHLHa17; Hpt; Scl; SCL/tal-1; tal-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011527
ORF Size:	987 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204834).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011527.2
RefSeq Size:	4161 bp
RefSeq ORF:	990 bp
Locus ID:	21349
UniProt ID:	P22091
Cytogenetics:	4 52.73 cM
Gene Summary:	Implicated in the genesis of hemopoietic malignancies. It may play an important role in hemopoietic differentiation. Serves as a positive regulator of erythroid differentiation. [UniProtKB/Swiss-Prot Function]



[View online »](#)