

Product datasheet for **MR204846**

Wdr42a (BC023704) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Wdr42a (BC023704) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Wdr42a
Synonyms:	H326, MGC38511
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR204846 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGGAGGACACGGGCCATTATTCATCAATGATGAAAGCCGAGGTCATGGCCACTCAGACGAAGAGGACG
 AGGAACAACCTCGGCACCGTGGACAGCGCAAACGGGCCAGCCGTGACCAGGACTCGTCAGATGATGAGCG
 AGCCCTGGAAGACTGGGTATCTTCAGAGACAACAGCCTTACCCGACCTCGCTGGCAAGCCCTTCCTGCG
 CTTGGGAGCGGGAAGTGGGTTCCAGTGCCCGCTTTGTGTACGAGGCTTGTGGAGCTAGAGTCTTTGTAC
 AGCGTTTCCGTTTGCAGCACGGGCTTGAGGGCCATACTGGTTGTGTCAATACCTTGCACTTTAACGAGCG
 TGGTACCTGGCTGGCCAGTGGCAGTGATGACCTGAAGGTGGTGGTGTGGGATTGGGTGCGGCGGCAGCCA
 GTACTGGACTTTGAGAGTGGCCACAAAAGTAATGTCTTCCAGGCCAAGTTCCTTCTTAACAGTGGTGACT
 CCACCCTGGCTATGTGTGCCCGGGATGGACAGGTCCGAGTCGAGAGCTGTCTGTACACAGTGTGCAA
 GAACACCAAGCGTGTGGCCAGCACAAGGGGGCATCTCACAAGTTGGCCCTAGAACCAGACTCCCCCTGC
 ACGTTCCTGTGAGCTGGTGAGGATGCCGTCTTACCATTGACCTGAGACAAGACCGGCCAGCTTCGA
 AACTAGTGGTGACAAAAGAGAAGGAGAAGAAAGTGGGGCTGTATACAATCTATGTGAACCCTGCCAATAC
 CCACCAAGTTGAGTGGGTGGAAGAGATCAGTATGTAAGGATTTATGACCAGAGGAAAATTGATGAGAAT
 GAGAACAACGGAGTCCTGAAAAAATCTGTCCTCACCATTGAGCTCCTGGCCAGTTACAATGATGAAGA
 CATTTACCTTTCAACTCTTCTCAGTGTATGGGGCCAGTATATTAAGAGATACAAAGGCCATAGAAAT
 AATGCTACAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)


Protein Sequence: >MR204846 protein sequence
 Red=Cloning site Green=Tags(s)

MEDTGHYSINDESIRGHGHSDEEDEEQPRHRGQRKRASRDQDSSDDERALEDWVSSETTALPRPRWQALPALRERELGSSARFVYEACGARVFVQRFRLQHGLEGHTGCVNTLHFNQRTWLASGSDDLKVVVWDWVRRQPVLD FESGHKSNVFQAKFLPNSGDSLAMCARDGQVRVAELSATQCCKNTRVAQHKGASHKLALPDSPCTFLSAGEDAVVFTIDLRQDRPASKLVVTKEKEKKVGLYTIYVNPANTHQF AVGGRDQYVRIYDQRKIDENENNGVLKKFCPHHLSSWPVTMMKTFSTLLTMGPSILRDTKAIEIMLQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC023704

ORF Size: 990 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC023704](#), [AAH23704](#)

RefSeq Size: 2276 bp

RefSeq ORF: 992 bp


Locus ID: 98193

Cytogenetics: 1 79.54 cM

MW: 37.4 kDa

Gene Summary: May function as a substrate receptor for CUL4-DDB1 E3 ubiquitin-protein ligase complex. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR204846