

Product datasheet for **MR204863L3V**

Alad (NM_008525) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Alad (NM_008525) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Alad
Synonyms:	ALADH; Lv
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008525
ORF Size:	993 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204863).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008525.3
RefSeq Size:	4536 bp
RefSeq ORF:	993 bp
Locus ID:	17025
UniProt ID:	P10518
Cytogenetics:	4 33.17 cM
Gene Summary:	Catalyzes an early step in the biosynthesis of tetrapyrroles. Binds two molecules of 5-aminolevulinate per subunit, each at a distinct site, and catalyzes their condensation to form porphobilinogen (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)