

## Product datasheet for MR204865

### Cd300lf (NM\_145634) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Cd300lf (NM_145634) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cd300lf
Synonyms:	CLIM1; CLM-1; CLM1; Digr2; F730004D16Rik; IgSF13; IREM1; LMIR3; Pigr3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR204865 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC

ATGCATTTGTCATTGCTGGTCCCCTTTCTCTTCTGGATCACAGGCTGCTGCACGGCTGAGGATCCAGTCA  
CAGGTCCAGAGGAGGTGAGCGGTCAGGAGCAGGGCTCCTTGACAGTGCAGTGCCGATACCTCAGGCTG  
GAAGGATTACAAGAAGTACTGGTGCCAAGGAGTTCCTCAGAGATCATGTAAGACTCTTGTTGAAACCGAT  
GCATCAGAGCAGCTGGTGAAGAAGAACCCTGTGCCATCAGGGACAACCAGAGAGACTTCATCTTCACAG  
TGACCATGGAGGATCTGAGGATGAGCGATGCTGGCATTACTGGTGTGGAATTACGAAAGGTGGACTTGA  
TCCCATGTTTAAAGTTACTGTGAACATTGGCCAGTGCCAACCATGCCCCCATCACCTCCACCACCACC  
ATCTTCACAGTGACAACCACAGTAAAAGAGACCAGCATGTTTCCAACGCTGACTAGCTACTCTGATA  
ACGGGCATGGCGGTGGTACAGTGGCGGTGGTGAAGATGGCGTCGGTGATGGGTTTCTGGATCTCAGTGT  
GCTCCTCCAGTCATCTCTGCAGTCCTGTTGCTTCTCCTGTTGGTGGCCTCGCTCTTTGCTTGGAGGATG  
GTGAGGAGACAGAAGAAAGCTGCTGGGCCACCATCAGAGCAGGCACAGTCTCTGGAGGGTGATCTCTGTT  
ATGCAGACCTGTCCCTGAAGCAGCCAGAACCTCCCCTGGCTCCTCTTGGAAAAAGGGCTCCTCCATGTC  
CTCCTCTGGCAAGGACCACCAAGAGGAAGTGAATATGTACCATGGCTCCCTTTCCAGGGAGGAGGTT  
TCATATGCCGCTCTGACTTTGGCCGGCTTGGGTCAGGAGCCTACTATGGCAATACTGGCTGCCCCATCA  
CCACGTTCCAGGACAGGCCTTGAAGAGGAGACCACAGAGTACAGCAGCATCAGGAGGCCCTTGCTGC  
AGCCATGCCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR204865 protein sequence  
Red=Cloning site Green=Tags(s)

MHLSLLVPFLFWITGCCTAEDPVTGPEEVSGQEQQSLTVQCRYTSGWKDYKKYWCQGVPPQRCKTLVETD  
 ASEQLVKKNRVSIIRDNRDFIFVTMEDLRMSDAGIYWCGITKGGLDPMFKVTVNIIGPVPTMPPITSTTT  
 IFVTVTTVKETSMPFTLTSYSDNGHGGGDSGGGEDGVGDGFLDLVLLPVISAVLLLLLLVASFVAVRM  
 VRRQKKAAGPPSEQAQSLQEGDLQCYADLSLKQPRTPSGSSWKKGSSMSSSGKDHQEEVEYVTPMAPFPREEV  
 SYAALTLAAGLQPEPTYGNTGCPITHVPRTGLEEETTEYSSIRRPLPAAMP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_145634

**ORF Size:** 993 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_145634.3](#)

**RefSeq Size:** 1864 bp

**RefSeq ORF:** 993 bp

**Locus ID:** 246746


**UniProt ID:** [Q6SIQ7](#)

**Cytogenetics:** 11 E2

**MW:** 36 kDa

**Gene Summary:** Acts as an inhibitory receptor for myeloid cells and mast cells (PubMed:17438331). Positively regulates the phagocytosis of apoptotic cells (efferocytosis) via phosphatidylserine (PS) recognition; recognizes and binds PS as a ligand which is expressed on the surface of apoptotic cells (PubMed:21865548). Plays an important role in the maintenance of immune homeostasis, by promoting macrophage-mediated efferocytosis and by inhibiting dendritic cell-mediated efferocytosis (PubMed:26768664). Negatively regulates Fc epsilon receptor-dependent mast cell activation and allergic responses via binding to ceramide which acts as a ligand (PubMed:23123064). May act as a coreceptor for interleukin 4 (IL-4). Associates with and regulates IL-4 receptor alpha-mediated responses by augmenting IL-4- and IL-13-induced signaling (PubMed:26124135). Negatively regulates the Toll-like receptor (TLR) signaling mediated by MYD88 and TRIF through activation of PTPN6/SHP-1 and PTPN11/SHP-2 (By similarity). Inhibits osteoclast formation (PubMed:14662855). Induces macrophage cell death upon engagement (PubMed:18097021).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR204865