

Product datasheet for **MR204868L4V**

Mcl1 (NM_008562) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mcl1 (NM_008562) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mcl1
Synonyms:	AW556805; Mcl-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008562
ORF Size:	993 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204868).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008562.3 , NP_032588.1
RefSeq Size:	3498 bp
RefSeq ORF:	996 bp
Locus ID:	17210
UniProt ID:	P97287
Cytogenetics:	3 40.74 cM
Gene Summary:	Involved in the regulation of apoptosis versus cell survival, and in the maintenance of viability but not of proliferation. Mediates its effects by interactions with a number of other regulators of apoptosis. Isoform 2 has antiapoptotic activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)