

## Product datasheet for **MR204870**

### Cdk10 (NM\_194444) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cdk10 (NM\_194444) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Cdk10  
**Synonyms:** BC017131  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >MR204870 representing NM\_194444  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

CTGGGAAGATGCCGAAGCGTTAAGGAGTTTGAAGAAGCTGAACCGGATTGGCGAGGGCACCTATGGCATCG  
TGTATCGGGCCAGGGATACCCAGACAGATGAAATTGTCGCCCTGAAGAAGGTGCGGATGGACAAAGAGAA  
GGATGGCATCCCCATCAGCAGCCTGCGTGAGATCACACTGCTTTCGCTCTCCGCCATCCAAACATTGTG  
GAGCTGAAGGAGGTGTTGTGGCAACACCTGGAGAGCATCTTCTGGTTCATGGTTACTGCGAAACAAG  
ATCTGGCCAGCCTATTGGAAAATATGCCAACACCTTCTCGGAGGCCAGGTTAAATGCATCATGCTACA  
GGTGTCTCGTGGCCTTACGTACCTGCACAGGAACCTTCATCATCCACAGGGACCTGAAGGTGTCCAATTG  
CTCATGACAGACAAGGGCTGTGTAAGACAGCTGATTTTGGCCTGGCTCGGGCCTATGGTGTCCAGTAA  
AGCCAATGACTCCCAAGTTGTTACCCTCTGGTACCGAGCCCCAGAGCTGCTGCTTGGAACTACCACCCA  
GACTACCAGCATTGACATGTGGGCTGTGCGGCTGCATCCTGGCAGAGCTGCTGGCCATAAGCCCCCTCTC  
CCTGGCACTTCCGAGATCCACCAGATCGACTTGATTGTACAGCTGTTGGGGACACCGAGTGAGAATATCT  
GGCCGGGTTTCTCCAAGCTGCCGCTGGCCGGCCAGTACAGCTTGAGGAAACAGCCCTATAACAACCTCAA  
GCACAAGTTCCTGGCTCTCAGAGGCCGACTCCGCTGCTCAACTTCTTTCATGTATGACCCCTAAG  
AAAAGGGCAACCTCGGGAGACTGCTGGAGAGCTCCTACTTCAAGGAGAAGCCCCCTGCCCTGTGAACCGG  
AGCTCATGCCTACCTTCCCCACCACCGCAATAAGCGTGTGCCAGCTGCCGCTGAAGGGCAGAGCAA  
ACGATGCCGGCC

**ACGCGT**ACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR204870 representing NM\_194444  
Red=Cloning site Green=Tags(s)

LGRCRSVKEFEKLNRIEGEGTYGIVYRARDTQTDEIVALKKVRMDKEKDGIPISLREITLLLRLRHPNIV  
 ELKEVVGNHLESIFLVMGYCEQDLASLLENMPTPFSEAQVKCIMLQVLRGLQYLHRNFIHRDLKVSNL  
 LMTDKGCVKTADFGLARAYGVPVKPMTPKVVTLYRAPELLGTTTQTTSIDMWAVGCILAELLAHKPLL  
 PGTSEIHQIDLIVQLLGTSPENIWPFGSKLPLAGQYSLRKQPYNNLKHKFPWLSEAGLRLLNLFMYDPK  
 KRATSGDCESSYFKEKPLPCEPELMPFPHHRNKRAAPAAAEGQSKRCRP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_194444

**ORF Size:** 993 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_194444.1</a>
<b>RefSeq Size:</b>	1623 bp
<b>RefSeq ORF:</b>	870 bp
<b>Locus ID:</b>	234854
<b>UniProt ID:</b>	<a href="#">Q3UMM4</a>
<b>Cytogenetics:</b>	8 E1
<b>MW:</b>	33.1 kDa
<b>Gene Summary:</b>	The protein encoded by this gene belongs to the CDK (cyclin-dependent kinase) subfamily of the Ser/Thr protein kinase family. The CDK subfamily members are highly similar to the gene products of <i>S. cerevisiae</i> cdc28, and <i>S. pombe</i> cdc2, and are known to be essential for cell cycle progression. The human ortholog has been shown to play a role in cellular proliferation. Multiple transcript variants encoding different isoforms have been found for this gene. A related pseudogene exists on chromosome 1. [provided by RefSeq, Jan 2010]