

Product datasheet for **MR204872L4V**

Lrrc26 (NM_146117) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lrrc26 (NM_146117) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lrrc26
Synonyms:	RP23-132N23.19
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_146117
ORF Size:	996 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204872).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_146117.2 , NP_666229.1
RefSeq Size:	1207 bp
RefSeq ORF:	996 bp
Locus ID:	227618
UniProt ID:	Q91W20
Cytogenetics:	2 A3


[View online »](#)

Gene Summary:

Auxiliary protein of the large-conductance, voltage and calcium-activated potassium channel (BK alpha). Required for the conversion of BK alpha channels from a high-voltage to a low-voltage activated channel type in non-excitabile cells. These are characterized by negative membrane voltages and constant low levels of calcium (By similarity).[UniProtKB/Swiss-Prot Function]