

Product datasheet for **MR204896L3V**

Npdc1 (NM_008721) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Npdc1
Synonyms:	AI314472; NPDC-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008721
ORF Size:	996 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR204896).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_008721.4 , NP_032747.2
RefSeq Size:	1522 bp
RefSeq ORF:	999 bp
Locus ID:	18146
UniProt ID:	Q64322
Cytogenetics:	2 A3



Gene Summary:

Suppresses oncogenic transformation in neural and non-neural cells and down-regulates neural cell proliferation. Might be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]