

## Product datasheet for **MR204908**

### Ifngr2 (NM\_008338) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Ifngr2 (NM_008338) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ifngr2
Synonyms:	Ifgr2; Ifgt
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR204908 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCGGCCTTTGCCACTGTGGCTGCCGTCGCTGCTCTGTGGGCTCGGTGCCGCGCGTCTCGCCAG  
ACTCGTTTTCCAGCTTGGCGCCCTCTGAACCAAGGCTTACCTGTACAATGATGAGCAGATTCTAAC  
TTGGGAGCCGTCACCTTCCAGCAATGACCCAAGACCAGTGGTCTACCAGGTGGAATATAGCTTCATCGAT  
GGCTCTTGGCATAGGTTGCTGGAGCCGAAGTACGGACATCACAGAGACAAAGTGTGACTTAACAGGAG  
GCGGCCGCTTGAAGCTTTCCACACCCATTACAGTCTTCTCGGGTGGCAGCAAGCGAGGAACT  
CACTTCCAAGTGGTGGGCTGGAGCCATTTCAACTATGAGAATGTTACTGTTGGACCTCCGAAAAAC  
ATCTCGGTGACCCAGGAAAAGTTCCCTCGTCATACACTTCTCCCTCCCTTGTATGTGTTCCACGGGG  
CACTTTTCAGTATCTTGTCCACTACTGGGAAAAGTCAGAAACCAACAGGAACAGGTTGAAGGCCCTTT  
CAAGAGCAACTCCATTGTGCTGGCAATCTGAAGCCATACAGAGTATATTGTTTACAAACTGAGGCACAA  
CTGATTTTGA AAAA AAAAAATCCGACCACATGGGCTCTTGAGCAATGTATCCTGTACGAAACACAG  
CAAATGCCTCCGCCAGGCTGCAGCAAGTCACTGATCCGTTGGGCATCTTCGATTGCTGCTCGGCCT  
GACGGGCGCCTGCTTACCCTGTTCTCAAATACCAAAGCCGAGTGAAGTACTGGTTTCAGGCTCCGCCA  
AACATCCCGAAACAAATCGAAGAGTATCTAAAGGACCCAGACCAATTCATCTTAGAGGCTTGGACAAGG  
ACGGTTCACCGAAGGAGGACTCCTGGGACTCCGTGTCAATTATTTCTTCCAGAAAAGGAGCGAGATGA  
TGTGCTCAAACACCG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)



**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_008338.4](#)

**RefSeq Size:** 2070 bp


**RefSeq ORF:** 999 bp

**Locus ID:** 15980

**Cytogenetics:** 16 53.07 cM

**MW:** 37.5 kDa

**Product images:**



Circular map for MR204908