

Product datasheet for **MR204913L3V**

Etfa (NM_145615) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Etfa (NM_145615) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Etfa |
| Synonyms: | 2010200I21Rik; D9Ertd394e |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_145615 |
| ORF Size: | 1002 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR204913). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_145615.2 |
| RefSeq Size: | 1464 bp |
| RefSeq ORF: | 1002 bp |
| Locus ID: | 110842 |
| UniProt ID: | Q99LC5 |
| Cytogenetics: | 9 B |



[View online »](#)

Gene Summary:

Heterodimeric electron transfer flavoprotein that accepts electrons from several mitochondrial dehydrogenases, including acyl-CoA dehydrogenases, glutaryl-CoA and sarcosine dehydrogenase. It transfers the electrons to the main mitochondrial respiratory chain via ETF-ubiquinone oxidoreductase (ETF dehydrogenase). Required for normal mitochondrial fatty acid oxidation and normal amino acid metabolism.[UniProtKB/Swiss-Prot Function]