

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for MR204934L4V

Gapdh (NM_008084) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gapdh (NM_008084) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gapdh
Synonyms:	Ga; Gapd
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008084
ORF Size:	999 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204934).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<u>NM 008084.2, NP 032110.1</u>
RefSeq Size:	1444 bp
RefSeq ORF:	1002 bp
Locus ID:	14433
UniProt ID:	<u>P16858</u>
Cytogenetics:	6 59.32 cM



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Gene Summary:This gene encodes a member of the glyceraldehyde-3-phosphate dehydrogenase protein
family. The encoded protein has been identified as a moonlighting protein based on its ability
to perform mechanistically distinct functions. The encoded protein was originally identified as
a key glycolytic enzyme that converts D-glyceraldehyde 3-phosphate (G3P) into 3-phospho-D-
glyceroyl phosphate. Subsequent studies have assigned a variety of additional functions to
the protein including nitrosylation of nuclear proteins, the regulation of mRNA stability, and
acting as a transferrin receptor on the cell surface of macrophage. Alternative splicing results
in multiple transcript variants. Many pseudogenes similar to this locus are found throughout
the mouse genome. [provided by RefSeq, Jan 2014]

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US