

Product datasheet for MR204941

Cops5 (NM_013715) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cops5 (NM_013715) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cops5
Synonyms:	AI303502; CSN5; Jab1; Mov34; Sgn5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR204941 representing NM_013715 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAGCTTCCGGGAGTGGTATGGCCAGAAAACCTGGGAATTGGCCAACAACATGCAGGAAGCGCAGA
GTATCGATGAAATCTACAAATATGACAAAAACAACAAGAAATCCTGGCGGCGAAACCTGGACTAA
GGATCACCCTACTTTAAATACTGCAAAATCTCAGCATTGGCTCTACTGAAAATGGTGATGCATGCCAGG
TCAGGAGGCAACTTGAAGTGATGGGTTTGTATGCTCGGAAAAGTCGACGGCGAGACCATGATCATCATGG
ACAGTTTCGCTTTCCTGTAGAGGGCACAGAACTCGAGTAAATGCTCAAGCTGCTGCGTATGAGTATAT
GGCTGCATACATAGAAAATGCCAAACAGGTTGGCCGCCTTGAGAAATGCAATCGGTTGGTATCATAGCCAC
CCTGGTTATGGCTGCTGGCTCTCCGGGATTGATGTTAGTACACAGATGCTGAACCAGCAGTTTCAAGAAC
CATTGTAGCAGTGATTGATCCAACCAACAATCTCTGCAGGAAAAGTGAATCTTGGCGCCTTTAG
GACATATCCAAAGGGCTACAAACCTCCTGATGAAGGACCTTCTGAGTACCAGACTATCCCCTTAATAAA
ATAGAAGATTTGGCGTGCCTGCAAAACAATATTATGCCTTAGAAGTCTCATATTTCAAATCATCTTTGG
ATCGTAAACTACTTGAGCTTTTGTGGAATAAATACTGGGTGAATACCCTGAGTTCCTCTAGCTTGCTTAC
TAATGCAGACTACACCACAGGCCAGGTGTTGATTTGTCTGAGAAGTTAGAGCAGTCGGAAGCCCAACTG
GGACGTGGCAGTTTCATGTTGGCTTAGAAACACATGACCGCAAGTCGGAAGACAACCTTGCCAAAGCTA
CTAGAGACAGCTGTAAAACCACCATAGAAGCCATCCATGGACTGATGTCTCAGGTTATTAAGGATAAACT
GTTAATCAGATTAACGTTGCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR204941 representing NM_013715
 Red=Cloning site Green=Tags(s)


MAASGSGMAQKTWELANNMQEAQSIDEIYKYDKKQQQEILAAKPWTKDHHYFKYCKISALALLKVMVHAR
 SGGNLEVMGLMLGKVDGETMIIMDSFALPVEGTETRVNAQAAAYEYMAAYIENAKQVGRLENAIGWYHSH
 PGYGCWLSGIDVSTQMLNQFQEPFVAVVIDPRTTISAGKVNLGAFRTYPKGYKPPDEGPSEYQTIPLNK
 IEDFGVHCKQYYALEVSYFKSSSLDRKLELLWNKYWVNTLSSSSLLTNADYTTGQVFDLSEKLEQSEAQL
 GRGSFMLGLETHDRKSEDKLAKATRDSCSTTIEAIHGLMSQVIKDKLFNQINVA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9034_e02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_013715


ORF Size: 1002 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_013715.2 , NP_038743.1
RefSeq Size:	1273 bp
RefSeq ORF:	1005 bp
Locus ID:	26754
UniProt ID:	O35864
Cytogenetics:	1 2.29 cM
MW:	38 kDa
Gene Summary:	Probable protease subunit of the COP9 signalosome complex (CSN), a complex involved in various cellular and developmental processes. The CSN complex is an essential regulator of the ubiquitin (Ubl) conjugation pathway by mediating the deneddylation of the cullin subunits of the SCF-type E3 ligase complexes, leading to decrease the Ubl ligase activity of SCF-type complexes such as SCF, CSA or DDB2. Promotes the proteasomal degradation of BRSK2. The complex is also involved in phosphorylation of p53/TP53, c-jun/JUN, I κ B α /NFKBIA, ITPK1 and IRF8, possibly via its association with CK2 and PKD kinases. CSN-dependent phosphorylation of TP53 and JUN promotes and protects degradation by the Ubl system, respectively. In the complex, it probably acts as the catalytic center that mediates the cleavage of Nedd8 from cullins. It however has no metalloprotease activity by itself and requires the other subunits of the CSN complex. Interacts directly with a large number of proteins that are regulated by the CSN complex, confirming a key role in the complex.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR204941