

Product datasheet for **MR204975L4V**

Lilrb4a (NM_013532) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lilrb4a (NM_013532) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lilrb4a
Synonyms:	CD85K; gp49; Gp49b; HM18; ILT3; Lilrb4; LIR-5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013532
ORF Size:	1008 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204975).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013532.2
RefSeq Size:	1497 bp
RefSeq ORF:	1008 bp
Locus ID:	14728
UniProt ID:	Q64281
Cytogenetics:	10 25.5 cM



[View online »](#)

Gene Summary:

Receptor for class I MHC antigens. Involved in the down-regulation of the immune response and the development of tolerance. Interferes with TNFRSF5-signaling and NF-kappa-B up-regulation. Inhibits receptor-mediated phosphorylation of cellular proteins and mobilization of intracellular calcium ions (By similarity).[UniProtKB/Swiss-Prot Function]