

Product datasheet for **MR204978L3V**

**Lenti ORF particles, MGC:31379 (Myc-DDK-tagged) - Mouse cDNA clone MGC:31379
 IMAGE:4240209, 200ul, >10⁷ TU/mL (Parp8) (BC022679) Mouse Tagged ORF Clone Lentiviral
 Particle**

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lenti ORF particles, MGC:31379 (Myc-DDK-tagged) - Mouse cDNA clone MGC:31379 IMAGE:4240209, 200ul, >10 ⁷ TU/mL (Parp8) (BC022679) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Parp8
Synonyms:	2810430O08Rik; ARTD16; D13Ertd275e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC022679
ORF Size:	1005 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204978).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC022679.1
RefSeq Size:	1536 bp
RefSeq ORF:	1007 bp
Locus ID:	52552
Cytogenetics:	13 65.32 cM


[View online »](#)

Gene Summary: Mono-ADP-ribosyltransferase that mediates mono-ADP-ribosylation of target proteins.
[UniProtKB/Swiss-Prot Function]