

Product datasheet for **MR204979L3V**

St6galnac5 (NM_012028) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	St6galnac5 (NM_012028) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	St6galnac5
Synonyms:	AI851940; Siat7e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012028
ORF Size:	1008 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204979).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012028.2
RefSeq Size:	2026 bp
RefSeq ORF:	1008 bp
Locus ID:	26938
UniProt ID:	Q9QYJ1
Cytogenetics:	3 H3
Gene Summary:	Involved in the biosynthesis of ganglioside GD1a from GM1b. It exhibits higher activity with glycolipids than with glycoproteins.[UniProtKB/Swiss-Prot Function]



[View online »](#)