

Product datasheet for **MR205011L4V**

Tmem59l (NM_182991) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmem59l (NM_182991) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tmem59l
Synonyms:	5330410G16Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_182991
ORF Size:	1014 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205011).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_182991.2
RefSeq Size:	1491 bp
RefSeq ORF:	1014 bp
Locus ID:	67937
UniProt ID:	Q7TNI2
Cytogenetics:	8 B3.3
Gene Summary:	Modulates the O-glycosylation and complex N-glycosylation steps occurring during the Golgi maturation of APP. Inhibits APP transport to the cell surface and further shedding (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)