

Product datasheet for **MR205104L4V**

Ankrd32 (BC009101) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ankrd32 (BC009101) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ankrd32
Synonyms:	2700017A04Rik; Ankrd32; AW545819; Brctd1; Brctx; C730024G01Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC009101
ORF Size:	1020 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205104).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC009101.1
RefSeq Size:	1601 bp
RefSeq ORF:	1022 bp
Locus ID:	105377
Cytogenetics:	13 C1



[View online »](#)

Gene Summary:

Plays a role in the DNA damage response (DDR) pathway by regulating postreplication repair of UV-damaged DNA and genomic stability maintenance. The SLF1-SLF2 complex acts to link RAD18 with the SMC5-SMC6 complex at replication-coupled interstrand cross-links (ICL) and DNA double-strand breaks (DSBs) sites on chromatin during DNA repair in response to stalled replication forks. Promotes the recruitment of SLF2 and the SMC5-SMC6 complex to DNA lesions.[UniProtKB/Swiss-Prot Function]