

Product datasheet for **MR205105L4V**

Lims2 (NM_144862) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Lims2
Synonyms:	PINCH2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_144862
ORF Size:	1023 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR205105).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_144862.1
RefSeq Size:	1688 bp
RefSeq ORF:	1026 bp
Locus ID:	225341
UniProt ID:	Q91XD2
Cytogenetics:	18 B1



Gene Summary:

Adapter protein in a cytoplasmic complex linking beta-integrins to the actin cytoskeleton, bridges the complex to cell surface receptor tyrosine kinases and growth factor receptors.
[UniProtKB/Swiss-Prot Function]