

Product datasheet for **MR205119L3V**

Qk (BC053426) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Qk
Synonyms:	Qkl, I500005P18
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC053426
ORF Size:	1023 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR205119).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	BC053426.1
RefSeq Size:	2208 bp
RefSeq ORF:	1025 bp
Locus ID:	19317
Cytogenetics:	17 7.75 cM



Gene Summary:

RNA-binding protein that plays a central role in myelination (PubMed:10864952, PubMed:11917126). Also required for visceral endoderm function and blood vessel development (PubMed:11892011, PubMed:16470614). Binds to the 5'-NACUAAY-N(1,20)-UAAAY-3' RNA core sequence (PubMed:16041388). Acts by regulating pre-mRNA splicing, mRNA export, mRNA stability and protein translation, as well as cellular processes including apoptosis, cell cycle, glial cell fate and development (PubMed:10535969, PubMed:12467586, PubMed:11297509, PubMed:11917126, PubMed:15568022). Required to protect and promote stability of mRNAs such as MBP and CDKN1B which promotes oligodendrocyte differentiation (PubMed:10535969, PubMed:15568022). Participates in mRNA transport by regulating the nuclear export of MBP mRNA (PubMed:12467586). May also play a role in smooth muscle development (PubMed:14706070).[UniProtKB/Swiss-Prot Function]