

## Product datasheet for **MR205123L3V**

### Jund (NM\_010592) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Jund (NM_010592) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Jund
Synonyms:	Jund1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010592
ORF Size:	1023 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205123).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_010592.4</a>
RefSeq Size:	2878 bp
RefSeq ORF:	1026 bp
Locus ID:	16478
UniProt ID:	<a href="#">P15066</a>
Cytogenetics:	8 34.15 cM



[View online »](#)

**Gene Summary:**

The protein encoded by this intronless gene is a member of the JUN family, and a functional component of the AP1 transcription factor complex. This protein has been proposed to protect cells from p53-dependent senescence and apoptosis. Alternative translation initiation site usage results in the production of different isoforms (PMID:12105216). [provided by RefSeq, Nov 2013]