

Product datasheet for MR205130

Hoxc10 (NM_010462) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hoxc10 (NM_010462) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Hoxc10
Synonyms:	A1506621; Hox-3.6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR205130 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACATGCCCTCGCAATGTAACCTCCGAACCTCGTACGCGGAGCCCTTGGCTGCGCCGGGAGGAGGAGAGC
GCTATAACCGTAACGCAGGAATGTATATGCAATCTGGGAGTGACTTCAACTGCGGGGTGATGAGGGGCTG
CGGGCTCGCGCCCTCTCTCCAAGAGGGACGAGGGAGGCAGCCCAAACCTAGCCCTCAACACCTACCCG
TCCTACCTCTCGCAGCTGGACTCCTGGGGCAGCCCAAAGCCGCCTACCGCCTGGAACAACCTGTTGGCA
GGCCTCTGTCCTCTGTTCTACCCACCTAGTGTCAAGGAGGAGAATGTCTGCTGCATGTACAGTGCAGA
GAAGCGGGCGAAAAGTGGCCCTGAGGCAGCTCTACTCCACCCCTGCCGGAGTCTTGCTTGGGGAG
CAGGAGGTACCTGTACCCAGCTACTACCGAGCCAGCCCGAGCTACTCCGCGCTGGACAAAACGCCCACT
GTGCTGGGGCCAACGAGTTTCAAGCCCCCTTTGAGCAGCGGGCCAGTCTCAACTCGCGCACCCGAACATCT
GGAATCGCCTCAGCTTGGGGCAAAGTGAGTTTTCTGAGACCCCAAGTCCGACAGCCAGACCCCACT
CCCAATGAGATCAAGACAGAGCAAAGCCTGGCGGGCCCAAAGCCAGCCCTCGGAGAGCGAGAAGGAAC
GGGCAAGACCGCAGACTCCAGTCCAGACACCTCGGATAACGAAGCTAAAGAGGAGATAAAGGCAGAAAA
CACCACAGGAAATTGGCTGACAGCAAAGAGCGGAAGGAAGAAGAGGTGCCCTATACTAAACACCAGACG
CTGGAATTGGAGAAAGAATTTCTGTTCAATATGATTTGACGCGAGAGCGCCCTGGAGATTAGCAAGA
CCATTAACCTTACAGACAGACAAGTCAAAATCTGGTTTTCAAATCGCAGAATGAAACTCAAGAAAATGAA
CCGAGAGAATCGGATCCGGAACTGACCTCCAATTTTAATTTACC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR205130 protein sequence
Red=Cloning site Green=Tags(s)

MTCP RNVTPNSYAEPLAAPGGGERYNRNAGMYMQSGSDFNCGVMRGCLAPSLSKRDEGGSPNLALNTYP
 SYLSQLDSWGDPKAAAYRLEQPVGRPLSSCSYPSPVKEENVCCMYSAEKRAKSGPEALYSHPLPESCLGE
 HEVPVPSYYRASPSYALDKTPHCAGANEFEAPFEQRASLNSRTEHLESPQLGGKVSFPEPTPKSDSQTPS
 PNEIKTEQSLAGPKASPSESEKERAKTADSSPDTSDNEAKEEIKAEENTTGNWLTAKSGRKKRCPYTKHQ
 TLELEKEFLFNMYLTRERRLEISKITINLTDRQVKIWFQNRMRKMKKMNRENRIRELT SNFNFT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_010462

ORF Size: 1029 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_010462.5](#)

RefSeq Size: 1914 bp

RefSeq ORF: 1029 bp

Locus ID: 209448


UniProt ID: [P31257](#)

Cytogenetics: 15 58.03 cM

MW: 38.2 kDa

Gene Summary: Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR205130