

Product datasheet for **MR205131L4V**

Cd2bp2 (NM_027353) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cd2bp2 (NM_027353) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cd2bp2
Synonyms:	1500011B02Rik; 2410024K20Rik; AW048183
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_027353
ORF Size:	1029 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205131).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_027353.2
RefSeq Size:	3361 bp
RefSeq ORF:	1029 bp
Locus ID:	70233
UniProt ID:	Q9CWK3
Cytogenetics:	7 F3
Gene Summary:	Involved in pre-mRNA splicing as component of the U5 snRNP complex that is involved in spliceosome assembly.[UniProtKB/Swiss-Prot Function]



[View online »](#)