

Product datasheet for **MR205135L3V**

Abhd4 (NM_134076) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Abhd4 (NM_134076) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Abhd4
Synonyms:	1110035H23Rik; Abh4; AI429574
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_134076
ORF Size:	1029 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205135).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_134076.1
RefSeq Size:	2406 bp
RefSeq ORF:	1068 bp
Locus ID:	105501
UniProt ID:	Q8VD66
Cytogenetics:	14 C2



[View online »](#)

Gene Summary:

Lysophospholipase selective for N-acyl phosphatidylethanolamine (NAPE). Contributes to the biosynthesis of N-acyl ethanolamines, including the endocannabinoid anandamide by hydrolyzing the sn-1 and sn-2 acyl chains from N-acyl phosphatidylethanolamine (NAPE) generating glycerophospho-N-acyl ethanolamine (GP-NAE), an intermediate for N-acyl ethanolamine biosynthesis (PubMed:16818490, PubMed:25853435). Hydrolyzes substrates bearing saturated, monounsaturated, polyunsaturated N-acyl chains (PubMed:16818490, PubMed:25853435). Shows no significant activity towards other lysophospholipids, including lysophosphatidylcholine, lysophosphatidylethanolamine and lysophosphatidylserine (PubMed:16818490).[UniProtKB/Swiss-Prot Function]