

Product datasheet for **MR20517L3V**

BC003331 (BC003331) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	BC003331
Synonyms:	MGC7150
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC003331
ORF Size:	1275 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR205171).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [BC003331](#), [AAH03331](#)

RefSeq Size: 1621 bp

RefSeq ORF: 1034 bp

Locus ID: 226499

Cytogenetics: 1 G1

Gene Summary: May play a role in the trafficking of a subset of G-protein coupled receptors.[UniProtKB/Swiss-Prot Function]



