

Product datasheet for **MR205172L4V**

Lpar6 (NM_175116) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lpar6 (NM_175116) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lpar6
Synonyms:	2610302I02Rik; P2ry5; P2y5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_175116
ORF Size:	1035 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205172).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_175116.2 , NP_780325.1
RefSeq Size:	2468 bp
RefSeq ORF:	1035 bp
Locus ID:	67168
UniProt ID:	Q8BMC0
Cytogenetics:	14 D3
Gene Summary:	Binds to oleoyl-L-alpha-lysophosphatidic acid (LPA). Intracellular cAMP is involved in the receptor activation. Important for the maintenance of hair growth and texture (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)