

Product datasheet for **MR205186L3V**

Apoh (NM_013475) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Apoh (NM_013475) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Apoh
Synonyms:	B2GPI; beta-2-GPI; beta2-GPI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013475
ORF Size:	1035 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205186).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013475.3
RefSeq Size:	1165 bp
RefSeq ORF:	1038 bp
Locus ID:	11818
UniProt ID:	Q01339
Cytogenetics:	11 71.8 cM
Gene Summary:	Binds to various kinds of negatively charged substances such as heparin, phospholipids, and dextran sulfate. May prevent activation of the intrinsic blood coagulation cascade by binding to phospholipids on the surface of damaged cells.[UniProtKB/Swiss-Prot Function]



[View online »](#)