

Product datasheet for **MR205211L4V**

Usf2 (NM_011680) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Usf2 (NM_011680) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Usf2
Synonyms:	bHLHb12; Usf-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011680
ORF Size:	1038 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205211).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011680.2 , NP_035810.1
RefSeq Size:	2419 bp
RefSeq ORF:	1041 bp
Locus ID:	22282
UniProt ID:	Q64705
Cytogenetics:	7 19.27 cM
Gene Summary:	Transcription factor that binds to a symmetrical DNA sequence (E-boxes) (5'-CACGTG-3') that is found in a variety of viral and cellular promoters.[UniProtKB/Swiss-Prot Function]



[View online »](#)