

Product datasheet for **MR205222L3V**

Ghitm (NM_078478) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ghitm (NM_078478) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ghitm
Synonyms:	1010001P14Rik; C77840; MICS1; PTD010; Tmbim5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_078478
ORF Size:	1041 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205222).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_078478.1 , NP_510963.1
RefSeq Size:	2906 bp
RefSeq ORF:	1041 bp
Locus ID:	66092
UniProt ID:	Q91VC9
Cytogenetics:	14 21.29 cM
Gene Summary:	Required for the mitochondrial tubular network and cristae organization. Involved in apoptotic release of cytochrome c (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)