

Product datasheet for MR205223

Mrpl38 (BC002319) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mrpl38 (BC002319) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Mrpl38
Synonyms:	MGC8310, MGC41062, MRP-L3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR205223 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCAACGAGGACATCGACGTGAGCAACCTGGAGCGGCTGGAGAAGTACCGGAGCTTCGAGCGGTACC
GCCCGGGCGGAGCAGGAGGCGCGCCCCGACTGGTGGCGGACCTACCGGGAGCATTTCTGTAAGGA
GACAGATCCCAAAGACAAAATTGACATTGGGCTACCCCCACCCAGAGTCAGTCGGACCCAAAAGCTACTG
GAGCGGAAGCACTTTCTGCGGAACTTCGTGCCAATGTGGAAGAAGAACGGGCAGCCCGCTCCGCACGG
CTAGCATCCCACTGGAGGCAGTCGGGGCAGAGTGGGAGAGGACGTGTGGGCCCTACCACAAGCAGCGTCT
GGCTGAGTACTATGGCCTCTACCGAGACCTGTTCCACGGTGCCACTTTTGTGCCCTGGGTCCCTCTACAC
GTGGCTTACGCCGTGGGAGAAGAGGACCTGATCCCGGTGTATCATGGCAATGAAGTGACTCCGACCGAGG
CTTCCCGGGCCCCAGAGGTGACCTATGAGGCAGACAAGGATTCCTTGTGGACTGCTATTCATCAACTT
GGACGGACACCTGCTGGAGCCGGACGCTGAGTATGTTCACTGGCTGCTAACCAACATCCCAAGCAACAGA
GTGGCTGAAGGACAGGAGACATGTCCTACCTCCCTCCTTTCCCTGCTCGAGGCTCTGGCTTCCACCGCT
TTGCCCTTCTGCTCTTCAAGCAGGACAAGCCAATCAACTTCTCAGAGGACACCCGTCCCTCACCTGCTA
CCAGCTGGCCCAGAGGACCTTCCGCACTTTTGATTTCTACAAGAGACACCAGGAAGCCATGACTCCAGCC
GGTCTGGCTTTTTCCAGTCCCGCTGGGATGACTCTGTCAACCCACCTTCCACCGCTTCTAGACATGC
GGGAACCTGTGTTTGAATTTGTACGGCCACCCCTTACCACCCAAACAGAAGCGCTTCCCCACGAGCA
GCCCCACGCTATCTGGACCGATACAGGGACAGTCATGAGCCACCTATGGAATCTAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR205223 protein sequence
Red=Cloning site Green=Tags(s)

MPNEDIDVSNLERLEKYRSFERYRRRAEQEARAPHWWRTYREHFVKETDPKDKIDIGLPPPRVSRTQKLL
 ERKHFLRELNRANVEEERAARLRTASIPLEAVRAEWERTCGPYHKQRLAEYYGLYRDLFHGATFVPWVPLH
 VAYAVGEEDLIPVYHGNEVTPTEASRAPEVTYEADKDSLWTLF INLDGHILLEPDAEYVHWLLTNIPSNR
 VAEGQETCPYLPFPARGSGFHRFAFLLFKQDKPINFSEDRPSPCYQLAQRTFRTFDFYKRHQEAMTPA
 GLAFFQCRWDDSVTHTFHQLDMREPVFVFRPPPYHPKQKRFPHQEPLRYLDRYRDSHEPTYGIY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: BC002319

ORF Size: 1038 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC002319](#), [AAH02319](#)

RefSeq Size: 1406 bp


RefSeq ORF: 1040 bp

Locus ID: 60441

Cytogenetics: 11 E2

MW: 41.3 kDa

Product images:



Circular map for MR205223