

Product datasheet for **MR205233L4V**

Elavl2 (NM_207686) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Elavl2 (NM_207686) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Elavl2 |
| Synonyms: | Hub; mel-N1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_207686 |
| ORF Size: | 1041 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR205233). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_207686.1 , NP_997569.1 |
| RefSeq Size: | 3741 bp |
| RefSeq ORF: | 1044 bp |
| Locus ID: | 15569 |
| Cytogenetics: | 4 42.36 cM |
| Gene Summary: | Binds RNA. Seems to recognize a GAAA motif. Can bind to its own 3' UTR, the FOS 3' UTR and the ID 3' UTR.[UniProtKB/Swiss-Prot Function] |



[View online »](#)