

Product datasheet for MR205235L3

Zfyve27 (NM_177319) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zfyve27 (NM_177319) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Zfyve27
Synonyms:	2210011N02Rik; 9530077C24Rik; AI426636; AI593546; AI835681
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205235).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_177319
ORF Size:	1245 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_177319.3 , NP_796293.2
RefSeq Size:	5636 bp
RefSeq ORF:	1248 bp
Locus ID:	319740
UniProt ID:	Q3TXX3
Cytogenetics:	19 C3
Gene Summary:	Key regulator of RAB11-dependent vesicular trafficking during neurite extension through polarized membrane transport (By similarity). Promotes axonal elongation and contributes to the establishment of neuronal cell polarity (PubMed:24251978). Involved in nerve growth factor-induced neurite formation in VAPA-dependent manner. Contributes to both the formation and stabilization of the tubular ER network. Involved in ER morphogenesis by regulating the sheet-to-tubule balance and possibly the density of tubule interconnections (By similarity). Acts as an adapter protein that facilitates the interaction of KIF5A with VAPA, VAPB, SURF4, RAB11A, RAB11B and RTN3 and the ZFYVE27-KIF5A complex contributes to the transport of these proteins in neurons. Can induce formation of neurite-like membrane protrusions in non-neuronal cells in a KIF5A/B-dependent manner (PubMed:21976701). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR205235L3