

Product datasheet for **MR205238L2V**

Cnr2 (NM_009924) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cnr2 (NM_009924) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cnr2
Synonyms:	CB-2; CB2; CB2-R
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_009924
ORF Size:	1044 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205238).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009924.2
RefSeq Size:	3865 bp
RefSeq ORF:	1044 bp
Locus ID:	12802
UniProt ID:	P47936
Cytogenetics:	4 D3
Gene Summary:	Heterotrimeric G protein-coupled receptor for endocannabinoid 2-arachidonoylglycerol mediating inhibition of adenylate cyclase. May function in inflammatory response, nociceptive transmission and bone homeostasis.[UniProtKB/Swiss-Prot Function]



[View online »](#)