

## Product datasheet for **MR205247L3V**

### **Cd72 (BC003824) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Cd72 (BC003824) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cd72
Synonyms:	CD72c; Ly-19; Ly-32; Ly-m19; Lyb-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC003824
ORF Size:	1044 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205247).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC003824</a> , <a href="#">AAH03824</a>
RefSeq Size:	1423 bp
RefSeq ORF:	1046 bp
Locus ID:	12517
Cytogenetics:	4 23.04 cM
Gene Summary:	Plays a role in B-cell proliferation and differentiation.[UniProtKB/Swiss-Prot Function]



[View online »](#)