

Product datasheet for **MR205262L4V**

BC004004 (NM_030561) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BC004004 (NM_030561) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	BC004004
Synonyms:	2400006G15Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_030561
ORF Size:	1047 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205262).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030561.2
RefSeq Size:	2809 bp
RefSeq ORF:	1047 bp
Locus ID:	80748
UniProt ID:	Q99KU6
Cytogenetics:	17 A3.3
Gene Summary:	Exhibits histone deacetylase (HDAC) enhancer properties. May play a role in cell cycle progression and wound repair of bronchial epithelial cells.[UniProtKB/Swiss-Prot Function]



[View online »](#)