

Product datasheet for **MR205267L3V**

Scrt1 (BC057562) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Scrt1 (BC057562) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Scrt1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC057562
ORF Size:	1044 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205267).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC057562 , AAH57562
RefSeq Size:	3498 bp
RefSeq ORF:	1046 bp
Locus ID:	170729
Cytogenetics:	15 D3
Gene Summary:	Transcriptional repressor that binds E-box motif CAGGTG (By similarity). Appears to function downstream of proneural bHLH proteins in promoting neural differentiation. [UniProtKB/Swiss-Prot Function]


[View online »](#)