

Product datasheet for **MR205282L4V**

Pdlim2 (NM_145978) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pdlim2 (NM_145978) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pdlim2
Synonyms:	4732462F18Rik; Slim
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145978
ORF Size:	1050 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205282).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145978.1
RefSeq Size:	1514 bp
RefSeq ORF:	1050 bp
Locus ID:	213019
UniProt ID:	Q8R1G6
Cytogenetics:	14 D2



[View online »](#)

Gene Summary:

Probable adapter protein located at the actin cytoskeleton that promotes cell attachment. Necessary for the migratory capacity of epithelial cells. Overexpression enhances cell adhesion to collagen and fibronectin and suppresses anchorage independent growth. May contribute to tumor cell migratory capacity (By similarity).[UniProtKB/Swiss-Prot Function]