

Product datasheet for **MR205289L4V**

Fbxo31 (NM_133765) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fbxo31 (NM_133765) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fbxo31
Synonyms:	1110003O08Rik; 2310046N15Rik; Fbx14; Fbxo14
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133765
ORF Size:	1050 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205289).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133765.2 , NP_598526.2
RefSeq Size:	4372 bp
RefSeq ORF:	1524 bp
Locus ID:	76454
UniProt ID:	Q3TQF0
Cytogenetics:	8 E1



[View online »](#)

Gene Summary:

Component of some SCF (SKP1-cullin-F-box) protein ligase complex that plays a central role in G1 arrest following DNA damage. Specifically recognizes phosphorylated cyclin-D1 (CCND1), promoting its ubiquitination and degradation by the proteasome, resulting in G1 arrest (By similarity).[UniProtKB/Swiss-Prot Function]