

Product datasheet for **MR205334**

Glcci1 (BC050802) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Glcci1 (BC050802) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Glcci1
Synonyms:	2310047L21Rik; A130036A18Rik; Fam117c; Gig18; Tssn1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR205334 representing BC050802 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGAGATAAGCCACTCAGACACCTAGTTGCTGGGCAGAAGAGGGAGCAGAAAAACGATCACATCAGC
 GCTCTGCATCATGGGAAGTGCTGATCAACTGAAAGAGCAGATTGCCAACTTAGGCAACAGCTACAACG
 CAGCAAGCAAAGTAGTCGGCATAGTAAAGAGAAGGACAGACAGTCACCTCTCCATGGCAACCACATAACA
 ATCAGTCATACTCAGGCTATTGGATCAAGATCAGTCCCTATGCCTCTGTCAAACATATCAGTGCCAAAT
 CATCTGTTTCCCGTGTACCCTGCAATGTAGAAGGAATAAGTCCTGAATTGAAAAGGTATTCATCAAAGA
 AAACAATGGAAAGGAAGAAGTGCCAAGCCATTGGATATACCAGATGGTCGAAGAGCTCCGCTCCCTGCT
 CACTACAGGAGCAGTAGTACTCGAAGCATAGATACCCAGACACCTTCTGTCCAAGAGCGCAGCAGTAGCT
 GCAGCAGCCACTCCCCTTGTGTGTCCCCATTTGTCTCCGGAATCCCAGGATGGAAGTCCTTGTCAAC
 AGAAGATTTGCTCTATGATCGTGATAAAGACAGTGGGAGTAGCTCACCGTTACCCAAGTATGCTTCATCT
 CCCAAACCAACAACAGCTACATGTTCAAACGGGAGCCCCAGAGGGATGTGAGCGGGTGAAGGTCTTTG
 AGGAAATGGCGTCTCGTCAGCCTATCTCGGCCCTCTCTTTTCATGTCTGACAAAAACAAGGTAAATTT
 CATCCCAACCGGATCAGCTTTCTGTCTGTAAACTTCTAGGCCCTCTCTTACCTGCCTCTGACCTGATG
 CTCAAGAACTCCCTAATTCTGGCCAGAGCTCGGCTCTGGCAACACTGACCGTAGAGCAGCTTTCTCTCC
 GGGTTTCTTTCACATCTCTCTGTATGACACCAGCACAGCCGACTCCCTGGAGCCCTCTGCCAGCAGCC
 ATCTCAGCAGCAGCAGCTCTACAGGACTTGCAGGTGGAGGAACAGCTCTCCACTCAGAACTATGTGATG
 ATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR205334 representing BC050802
 Red=Cloning site Green=Tags(s)

MRDKATQTPSCWAEEGAERSHQRSASWGSADQLKEQIAKL RQQLQRSKQSSRHSKEKDRQSPLHGNHIT
 ISHTQAIGSRVPMPLSNISVPKSSVSRVPCNVEGISPELEKVF IKENNGKEEVSKPLDIPDGRRAPLPA
 HYRSSSTRSIDTQTPSVQERSSSCSSHSPCVSPFCPPESQD GSPCSTEDLLYDRDKDSGSSSPLPKYASS
 PKPNNSYMFKREPPEGCERVKVFEEMASRQPI SAPLFSCPDKNKVNFIPTGSAFCPVKLLGPLLPASDLM
 LKNSPNSGQSSALATLTVEQLSSRVSFTLSDDTSTADSLEPSAQ QPSQQQQLQLQVEEHVSTQNYVM
 I



TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: BC050802

ORF Size: 1053 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC050802.1](#)

RefSeq Size: 1782 bp


RefSeq ORF: 1055 bp

Locus ID: 170772

Cytogenetics: 6 A1

MW: 38.5 kDa

Product images:



Circular map for MR205334