

Product datasheet for **MR205388L3V**

Pon2 (NM_183308) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pon2 (NM_183308) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pon2
Synonyms:	6330405I24Rik; A1481612
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_183308
ORF Size:	1065 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205388).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_183308.2
RefSeq Size:	1850 bp
RefSeq ORF:	1065 bp
Locus ID:	330260
UniProt ID:	Q62086
Cytogenetics:	6 1.99 cM
Gene Summary:	Capable of hydrolyzing lactones and a number of aromatic carboxylic acid esters. [UniProtKB/Swiss-Prot Function]



[View online »](#)