

Product datasheet for **MR205418L4V**

Bag1 (NM_001171739) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bag1 (NM_001171739) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bag1
Synonyms:	BAG-1; Ra; Rap46
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001171739
ORF Size:	1065 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205418).
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001171739.1
RefSeq Size:	1365 bp
RefSeq ORF:	660 bp



[View online »](#)

Locus ID: 12017

UniProt ID: [Q60739](#)

Cytogenetics: 4 A5

Gene Summary: The oncogene Bcl2 encodes a membrane protein that blocks a step in a pathway leading to apoptosis or programmed cell death. The protein encoded by this gene binds to Bcl2 protein and is referred to as Bcl2-associated athanogene. It enhances the anti-apoptotic effects of Bcl2 and represents a link between growth factor receptors and anti-apoptotic mechanisms. At least two protein isoforms are encoded by this mRNA through the use of a non-AUG (CUG) start site and an alternative, downstream, AUG translation initiation site. [provided by RefSeq, Jul 2008]