

Product datasheet for MR205431

Oas1e (NM_145210) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Oas1e (NM_145210) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Oas1e
Synonyms:	C87484; oasl7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR205431 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC

ATGGCGAGGGAACCTTTCAGAACCCCAATCTGGAGGCTGGATAACTTCATAGAGGATCAACTCCTTCCTGACACCACCTTCTTACTGAGCTCAGAGCAGACATCGACTCCATAAGTGCTTTCCTGAAGGAGAGATGCTTCCAAGGTGCCACCCACCCTATGAGGGTCGCCAGGGTTGTGATGGGAGGCTCCTATGATGAACACAATGCACTCAAGGGAAGGTCAGAGGCCAACATGGTGGTGTCTTTAATGATCTCACCAGCTTTGAGGACCAGTTAAAGTGGCAGCAAGTGTTCATTGAAGAAATTCGGAACACCTGCTCCAGTTGCAGCAAGAGAAGCCATGTAAACTCAAGTTTGGAGTGCAGAGCTCAGAGGAGCCCAACACCAGGTCTCTGACCTTCAAGCTGTGCTCCCCAGAGCGCCAGCAGGAGGTGGAATTTGATGTGCAGCCAGCCTATGATGCTCTGTATGAAGGGGGATACTGCAAGTCTTTGAATCCATTAACACAAAGTCTACGCCAACTCATCCATGAGTGCACCACCCTGGAGAAAGGAGGGCGAGTTCCTATCTGCTTACCAGCCTTCATCAGAGCTTCTGAGGTATCGTGCGCCCAAGCTCTGGAACCTCATCCGTCTGGTCAAGCACTGGTATCAACTGTGTAAGGAGAAGCTGAGGGGACCCTGCCTCACAGTATGCCCTGGAGCTGCTCACAGTCTACGTCTGGGAATTTGGGGTCCAAGACAGCTTTGGACTCCATGCAGCCAGTGTCCGAACGGTCTTAGAACTGGTACCAAGTACAAAATGCCTTCTAATCTACTGGACATGGTTTTATGACTTTCGACCTGAGATCTCTGACTACCTGCAGGACAGATCAAAAAGCCAGGCCTCTGATCCTGGATCCGGCAGACCAACAAGGAACGTGGCTGGTTACAGACTTACAGGCATGGGACCTGCTGGCAAAAGGAGCTCAGACCTGGATACATTCAAATTTTTTTAGGAACTGTGATATGTCCTTGTGAATGGCTGGGAAAGTGTGCCAGAGAAACA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR205431 protein sequence
 Red=Cloning site Green=Tags(s)

MARELFRTPIWRLDNFIEDQLLPDPTTFL TELRADIDSISAF LKERCFQGATHPMRVARVVMGGSYDEHNA
 LKGRSEANMVVFFNDLTSFEDQLKWQV FIEEIRKHL LQLQEKPKLKFEVQSSEEPNTRSLTFKLCSP
 ERQQEVEFDVQPAYDAL YEGGYCKSFESINYNKVYAQLIHECTTLEKEGEFSICFDLHQSF LRYRAPKL
 WNLIRLVKHVYQLCKEKL RGPLPQYALELLTVVWVEFGVQDSFGLHAAQCFRTVLELVTKYKCLLIYWT
 WFYDFRPEISDYLHGQIKKARPLILD PADPTRNVAGSDLQAWDLLAKEAQTWIHSNFFRNCDMSLVNGWE
 VSPEKQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_145210

ORF Size: 1071 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_145210.2

RefSeq Size: 1715 bp


RefSeq ORF: 1071 bp

Locus ID: 231699

Cytogenetics: 5 60.64 cM

MW: 42 kDa

Product images:



Circular map for MR205431