

Product datasheet for **MR205439L4V**

Fam134b (Retreg1) (NM_025459) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Fam134b (Retreg1) (NM_025459) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Retreg1 |
| Synonyms: | 1810015C04Rik; AU015349; Fam134b; Retreg1d |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_025459 |
| ORF Size: | 1071 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR205439). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_025459.2 |
| RefSeq Size: | 2929 bp |
| RefSeq ORF: | 1071 bp |
| Locus ID: | 66270 |
| UniProt ID: | Q8VE91 |
| Cytogenetics: | 15 B1 |



[View online »](#)

Gene Summary:

Endoplasmic reticulum-anchored autophagy receptor that mediates ER delivery into lysosomes through sequestration into autophagosomes (PubMed:26040720). Promotes membrane remodeling and ER scission via its membrane bending capacity and targets the fragments into autophagosomes via interaction with ATG8 family proteins (PubMed:26040720). Required for long-term survival of nociceptive and autonomic ganglion neurons (PubMed:19838196, PubMed:26040720).[UniProtKB/Swiss-Prot Function]