

Product datasheet for **MR205452L3V**

Sord (NM_146126) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Sord (NM_146126) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Sord |
| Synonyms: | Sdh-1; Sdh1; Sodh-1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_146126 |
| ORF Size: | 1074 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR205452). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_146126.1 |
| RefSeq Size: | 2259 bp |
| RefSeq ORF: | 1074 bp |
| Locus ID: | 20322 |
| UniProt ID: | Q64442 |
| Cytogenetics: | 2 60.59 cM |



[View online »](#)

Gene Summary:

Polyol dehydrogenase that catalyzes the reversible NAD(+)-dependent oxidation of various sugar alcohols (By similarity). Is active with D-sorbitol (D-glucitol) leading to the C2-oxidized product D-fructose (PubMed:6852349). Is a key enzyme in the polyol pathway that interconverts glucose and fructose via sorbitol, which constitutes an important alternate route for glucose metabolism (By similarity). May play a role in sperm motility by using sorbitol as an alternative energy source for sperm motility and protein tyrosine phosphorylation (PubMed:18799757). Has no activity on ethanol. Cannot use NADP(+) as the electron acceptor (PubMed:6852349).[UniProtKB/Swiss-Prot Function]