

Product datasheet for MR205458L3

G6pc (NM_008061) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	G6pc (NM_008061) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	G6pc
Synonyms:	AW107337; G6P; G6Pase; G6pc1; G6pt; Glc-6-; Glc-6-Pase
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205458).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_008061
ORF Size:	1074 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_008061.3
RefSeq Size:	2245 bp
RefSeq ORF:	1074 bp
Locus ID:	14377
UniProt ID:	P35576
Cytogenetics:	11 D
Gene Summary:	<p>The enzyme encoded by this gene is a multisubunit integral membrane protein of the endoplasmic reticulum that is composed of a catalytic subunit and transporters for glucose-6-phosphate, inorganic phosphate, and glucose. This gene is one of three glucose-6-phosphatase catalytic-subunit-encoding genes in mouse. Glucose-6-phosphatase catalyzes the hydrolysis of D-glucose 6-phosphate to D-glucose and orthophosphate and is a key enzyme in glucose homeostasis, functioning in gluconeogenesis and glycogenolysis. Mutations in this gene cause glycogen storage disease type I (GSD1). [provided by RefSeq, Sep 2015]</p>

Product images:



Circular map for MR205458L3