

## Product datasheet for **MR205471L3V**

### Chad (NM\_007689) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Chad (NM_007689) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Chad
Synonyms:	SLRR4A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007689
ORF Size:	1077 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205471).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_007689.2</a>
RefSeq Size:	1653 bp
RefSeq ORF:	1077 bp
Locus ID:	12643
UniProt ID:	<a href="#">O55226</a>
Cytogenetics:	11 58.9 cM



[View online »](#)

**Gene Summary:**

Promotes attachment of chondrocytes, fibroblasts, and osteoblasts. This binding is mediated (at least for chondrocytes and fibroblasts) by the integrin alpha(2)beta(1). May play an important role in the regulation of chondrocyte growth and proliferation (By similarity).  
[UniProtKB/Swiss-Prot Function]