

Product datasheet for MR205485L4

Cxcr4 (NM_009911) Mouse Tagged Lenti ORF Clone

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Expression Plasmids
Product Name:	Cxcr4 (NM_009911) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Cxcr4
Synonyms:	b2b220Clo; CD184; Cmkar4; LESTR; PB-CKR; PBSF/SDF-1; Sdf1r
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205485).
Restriction Sites:	Sgfl-Mlul
Cloning Scheme:	
	Cloning sites used for ORF Shuttling:
	Sgf I ORF Mlu I GCG ATC GCC ATG // NNÑ ACG CGT



ACCN: ORF Size: NM_009911 1077 bp



View online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	xcr4 (NM_009911) Mouse Tagged Lenti ORF Clone – MR205485L4
OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <u>custsupport@origene.com</u> or by calling 301.340.3188 option 3 for pricing and delivery.
	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Met	 thod: 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM 009911.3</u>
RefSeq Size:	1821 bp
RefSeq ORF:	1080 bp
Locus ID:	12767
UniProt ID:	<u>P70658</u>
Cytogenetics:	1 56.43 cM

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

GRIGENE Cxcr4 (NM_009911) Mouse Tagged Lenti ORF Clone – MR205485L4

Gene Summary:Receptor for the C-X-C chemokine CXCL12/SDF-1 that transduces a signal by increasing
intracellular calcium ion levels and enhancing MAPK1/MAPK3 activation (PubMed:8962122,
PubMed:9295051, PubMed:9103415). Involved in the AKT signaling cascade (By similarity).
Plays a role in regulation of cell migration, e.g. during wound healing. Acts as a receptor for
extracellular ubiquitin; leading to enhanced intracellular calcium ions and reduced cellular
cAMP levels. Binds bacterial lipopolysaccharide (LPS) et mediates LPS-induced inflammatory
response, including TNF secretion by monocytes (By similarity). Involved in hematopoiesis and
in cardiac ventricular septum formation (PubMed:9634237, PubMed:9634238,
PubMed:9689100). Also plays an essential role in vascularization of the gastrointestinal tract,
probably by regulating vascular branching and/or remodeling processes in endothelial cells
(PubMed:9634237). Involved in cerebellar development. In the CNS, could mediate
hippocampal-neuron survival (PubMed:9634238, PubMed:9689100).[UniProtKB/Swiss-Prot
Function]

Product images:



Circular map for MR205485L4

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US





Double digestion of MR205485L4 using Sgfl and Mlul

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US