

Product datasheet for **MR205531L3V**

Adk (NM_134079) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adk (NM_134079) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Adk
Synonyms:	2310026J05Rik; 5033405D03Rik; AI255373; AI987814; Ak
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_134079
ORF Size:	1083 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205531).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_134079.1
RefSeq Size:	1737 bp
RefSeq ORF:	1086 bp
Locus ID:	11534
UniProt ID:	P55264
Cytogenetics:	14 A3
Gene Summary:	ATP dependent phosphorylation of adenosine and other related nucleoside analogs to monophosphate derivatives. Serves as a potential regulator of concentrations of extracellular adenosine and intracellular adenine nucleotides.[UniProtKB/Swiss-Prot Function]



[View online »](#)