

Product datasheet for MR205539

Ybx3 (NM_139117) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ybx3 (NM_139117) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ybx3
Synonyms: Csda; dbpA; Dpba; MSY3; MSY4; oxyR; Yb2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR205539 representing NM_139117
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGCGAGGGCGGGCAGGGCCACCACGGCGGCACCACGCTCCCGCAGGCCGGCGCCGACGCGCCCGCG
 CGGCGCCCCGGACCCCGCGCCTAAGAGCCCGCGGCCAGCGCGCGCCACGGCCCCGGCGCCCGCCG
 GCTGCTCGCGGGAGCCCGCGGAGACGCAGCCCCGGGCCCGCCCGGCTCATCAGCCCCCGGGG
 GCGGAGGACGCGGAGAAGAAAGTTCTCGCCACCAAGTCTTGGCACTGTCAAATGGTTCAACGTCAGAA
 ATGGATATGGATTTATAAACCGAAATGACACCAAGAAGATGTGTTGTACACCAGACGGCCATCAAGAA
 AAACAACCCACGCAAGTATCTGCGCAGTGTGGGGATGGAGAACTGTAGAGTTTGATGTAGTTGAAGGA
 GAAAAGGGTGCTGAAGCAGCAAACGTGACTGGCCAGATGGAGTTCTGTAGAAGGGAGTCGCTATGCTG
 CTGACCGGCGCCGGTACAGACGCGGCTACTATGGCAGGCGCCGTGGACCTCCCCGTAATTACGCTGGGGA
 GGAGGAGGAGGAAGGGAGCGGCAGCAGTGAAGGATTTGAGCCCCCTGCCGAGATGGGCAGTTCTCTGGG
 GCCAGGAATCAGCTGCGCCGCCCCAGTATCGCCCTCCATACCGGCAGCGGCGTTTCCCGCCTTACCACG
 TGGGACAGACCTTTGACCGTCGCTCACGGGTCTTCCCATCCCAACAGAATGCAGGCTGGTGAGATTGG
 GGAGATGAAGGATGGAGTCCCCGAGGGAACACAGCTCCAGGCTCATCGGAATCCGACTTACCGCCCGAGG
 TTCCGCGAGGGACCTGCTCGCCACGACCTGCCCTGCTATTGGAGAGGCTGAAGACAAAGAAAACCAGC
 AAGCGGCTAATGGTCAAACAGCCGTCTGCCCGCGTGGATTCCGACGTCCTACAACATAGGGCGCCG
 TTCCCGTCCCCTCAACGCTGTTTCAAGATGGCAAAGAGACCAAGGCAGGTGAAGCACCAACTGAGAAC
 CCCGCTCCAGCTACCGAACAGAGCAGTCCGAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR205539 representing NM_139117
Red=Cloning site Green=Tags(s)

MSEAGEATTGGTTLPQAAADAPAAAPPDPAPKSPAASGAPQAPAPAALLAGSPGGDAAPGPAPASSAPAG
 GEDAEEKVLATKVLGTVKWFNVRNGYGF INRNDTKEDVFVHQTAIKKNPRKYLRSVGDGETVEFDVVEG
 EKGAEANVTGPDGVPVEGSRYAADRRRYRRGYGRRRGPPRNYAGEEEEEEGSGSSEGFEPAAADGGQFSG
 ARNQLRRPQYRPPYRQRFFPPYHVGQTFDRRSRVFPHPNRMQAGEIGEMKDGVPETQLQAHRNPTYRPR
 FRRGPARPRPAPAIGEAEDKENQQAANGPNQPSARRGFRFPYNYRRRSRPLNAVSQDGKETKAGEAPTEN
 PAPATEQSSAE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja1819_g03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_139117


ORF Size: 1083 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

- OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
- RefSeq:** [NM_139117.2](#), [NP_620817.2](#)
- RefSeq Size:** 1862 bp
- RefSeq ORF:** 1086 bp
- Locus ID:** 56449
- UniProt ID:** [Q9JKB3](#)
- Cytogenetics:** 6 F3
- MW:** 39.3 kDa
- Gene Summary:** Binds to the GM-CSF promoter. Seems to act as a repressor (By similarity). Binds also to full-length mRNA and to short RNA sequences containing the consensus site 5'-UCCAUCA-3'. May have a role in translation repression.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR205539