

## Product datasheet for **MR205539L3V**

### Ybx3 (NM\_139117) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ybx3 (NM_139117) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ybx3
Synonyms:	Csda; dbpA; Dpba; MSY3; MSY4; oxyR; Yb2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_139117
ORF Size:	1083 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205539).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_139117.2</a> , <a href="#">NP_620817.2</a>
RefSeq Size:	1862 bp
RefSeq ORF:	1086 bp
Locus ID:	56449
UniProt ID:	<a href="#">Q9JKB3</a>
Cytogenetics:	6 F3
Gene Summary:	Binds to the GM-CSF promoter. Seems to act as a repressor (By similarity). Binds also to full-length mRNA and to short RNA sequences containing the consensus site 5'-UCCAUCA-3'. May have a role in translation repression.[UniProtKB/Swiss-Prot Function]



[View online »](#)