

Product datasheet for **MR205581L3V**

Pgbd1 (BC089360) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pgbd1 (BC089360) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pgbd1
Synonyms:	4921509E05Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC089360
ORF Size:	1089 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205581).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC089360 , AAH89360
RefSeq Size:	2980 bp
RefSeq ORF:	1091 bp
Locus ID:	319207
Cytogenetics:	13 A3.1



[View online »](#)